

BNP-32 ,porcine

PRODUCT DESCRIPTION

BNP-32 ,porcine is a 53-amino-acid synthetic peptide. A synthetic peptide used as a research tool in biochemistry, pharmacology, and biomedical research. This peptide enables the study of protein-protein interactions, receptor binding, enzyme substrate specificity, and signal transduction pathways. Available in high purity (>=95% by HPLC) as a lyophilized powder for reconstitution in appropriate buffers.

PRODUCT SPECIFICATIONS

Catalog Number	BPT2151
CAS Number	117345-87-6
Molecular Weight	3570.0863 Da
Molecular Formula	C149H250N52O44S3
Purity	Peak Area by HPLC >=95%
Form	Lyophilized
Source	Synthetic (Chemical Synthesis, SPPS)
Endotoxin	N/A (Synthetic peptide)

AMINO ACID SEQUENCE

One-Letter Code:

SPKTMRDSGCFGRRRLDRIGSLSGLCNVLRRY(Disulfidebridge:Cys10-Cys26)

Three-Letter Code:

H-Ser-Pro-Lys-Thr-Met-Arg-Asp-Ser-Gly-Cys-Phe-Gly-Arg-Arg-Leu-Asp-Arg-Ile-Gly-Ser-Leu-Ser-Gly-Leu-Gly-Cys-Asn-Val-Leu-Arg-Arg-Tyr-OH(Disulfidebridge:Cys10-Cys26)

RESEARCH APPLICATIONS

Application	Receptor binding
Biomarker Target	Natriuretic Peptide Receptors
Research Area	Cardiovascular
Sub-category	Natriuretic Peptides & Heart Failure

STORAGE & HANDLING

Storage Temperature	-20 deg C
Form	Lyophilized
Stability	Lyophilized product stable for 24 months at -20 deg C
Reconstitution	Reconstitute in sterile water or appropriate buffer
After Reconstitution	Aliquot and store at -20 to -80 deg C
Avoid	Repeated freeze-thaw cycles

AVAILABLE SIZES

Size Options	1 mg, 10 mg
--------------	-------------

Purity Method	HPLC, Mass Spectrometry (MS)
Identity	Confirmed by ESI-MS or MALDI-TOF
Purity	Peak Area by HPLC >=95%
Documentation	Certificate of Analysis (CoA) available upon request